Ra Iv K		
Longitud (mm)	Diâmetro (mm)	Número de catálogo
500	6	022-0062
	7	022-0072
	8	022-0082
	9	022-0092
	10	022-0102
	12	022-0122
	14	022-0142
	16	022-0162
	18	022-0182
	20	022-0202
	22	022-0222
	24	022-0242
	25	022-0252
	26	022-0262
	28	022-0282
	30	022-0302
	32	022-0322

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INSTRUCTIONS FOR USE

Knitted vascular prosthesis with collagen - Ra K

PRODUCT DESCRIPTION

The Ra K knitted vascular prosthesis is a tube made of polyester filament yarn using a warp-knitting technology. The front wall of the prosthesis is provided with a colour guiding line that makes it easier for the surgeon to orientate himself in the surgical wound during implantation. A continuous layer of chemically modified bovine collagen is applied to the surface of the prosthesis. It facilitates ingrowth of tissue into the implant and impermeability of the prosthesis wall at the time of implantation. Ra 1v K is straight, crimped, with collagen, Ra bv K is bifurcated, crimped, with collagen. Thanks to the crimping the graft retains a lumen of circular shape even when it is bent.

TYPE

Ra 1v K – knitted vascular prosthesis, straight, crimped, with collagen
Ra bv K – knitted vascular prosthesis, bifurcated, crimped, with collagen

INTENDED USE

Knitted vascular prostheses impregnated with bovine collagen are suitable for arterial reconstructions where application of an artificial vascular graft of corresponding lumen is indicated. These types of prostheses can be used in aorto-femoral, ilio-femoral regions, femoro-popliteal region above the knee joint, in the region of extra-cranial arteries and for extraanatomical bypasses. Using these prostheses, arterial reconstructions can be carried out without great loss of blood.

KONTRAINDICATIONS

- with persons provably allergic to bovine proteins
- in case of an increased hazard of infection
- for pediatric patients, can be used only in urgent cases, after individual consideration

NOTES AND PREPARATION OF THE PROSTHESIS

FOR IMPLANTATION

- The prosthesis is used directly without the necessity of any preparation before implantation.
- The prosthesis can be adjusted with scissors at the ends; it is strong, elastic, easy to pierce, and it does not fray.
- The prosthesis is implanted under full heparinisation of the patient.
- After the implantation, cancel the effect of heparinisation by means of neutralisation by administration of an appropriate dose of Protamine Sulfate i.v. (intravenous).
- The prosthesis must not be injured with sharp instruments in order not to damage its impregnation.
- If it is necessary to stop the flow in the prosthesis use „shod clamps“, i.e. clamps with a rubber tube pulled on the branches of the clamp.
- When handling the prosthesis ensure careful handling, protect from impact in order not to brake the integrity of impregnation.

IMPLANTATION TECHNIQUE

- The prosthesis is to be implanted under moderate tension, when performing anastomoses use nonabsorbable monofilament atraumatic sutures designed for vascular surgery.

- We recommend gradual intake of blood flow into the prosthesis.
- After the upper (central) anastomosis has been performed, release the proximal clamp on the artery for one or two pulse waves. Then close the clamp on the artery again and leave the blood column in the graft. If the upper anastomosis does not bleed it is possible to apply a shod clamp on the graft right below the anastomosis. If the prosthesis becomes covered with blood-dew or if it leaks blood, then, after a short while, once more let one or two pulse waves into the graft so
- that a blood column may be formed in the graft. This having done, it is possible to draw the graft through the tunnels to the distal anastomosis.
- Before you perform the lower (peripheral) anastomosis, make sure of the right tension of the graft and check whether the graft is not kinked or twisted somewhere. Before you perform the distal anastomosis, aspirate the blood gently from the graft using a suction apparatus, and perform the distal anastomosis. Before final knotting of the suture, leave a free space in the anastomosis and release the proximal clamp for a short while so that the prosthesis may be deaerated or, as the case may be, the thrombi may be flushed away. Then close again the proximal clamp on the prosthesis, tighten the sutures, and knot them. Then it is possible to remove gradually the clamps from the periphery and, finally, from the upper anastomosis area.
- In case of slight bleeding from the anastomoses area or of blood leakage through the wall of the graft it is possible to apply towels, moistened with physiological saline whose temperature must not exceed the body temperature, to the graft (for a period of 5 - 10 min.) while the heparine is neutralized at the same time.

STERILIZATION

The product is delivered sterilized by irradiation.
The product is sterile provided that its package has not been damaged or opened.

The product is intended for one use only.

THE PRODUCT CANNOT BE RE-STERILIZED BY NO MEANS!!

NOTE:

COUNCIL DIRECTIVE 93/42/EEC, Annex No. I, article 8.3., medical devices delivered in a sterile state must be designed, manufactured and packed in a non-reusable pack and/or according to appropriate procedures to ensure that they are sterile when placed on the market and remain sterile, under the storage and transport conditions laid down, until the protective packaging is opened or damaged. For this reason the medical device must not be reused because opening of the original pack (nonreusable pack) means damage to the manufacturer's original pack. Once the original pack has been opened and the product taken out of this pack the sterility of the product is broken. Once the original pack has been opened the product must not be repackaged and re-sterilized. The manufacturer cannot guarantee that a contamination of the product by microorganisms will not occur after removal of the product from its original pack and during handling.

EXPIRATION PERIOD

Three years from the date of sterilization. The expiration date is given on the package.

STORAGE AND HANDLING CONDITIONS

The products shall be stored in original packing in rooms designed for storing. The temperature in storage rooms must be 15 °C to 25 °C. The products must be protected against effects of direct sunlight. Vascular prosthesis with collagen yellows which is an accompanying characteristic of collagen. The yellowing does not affect the performance properties of the vascular prostheses with collagen. The vascular prosthesis must not freeze because it would cause crackling of the collagen layer. At an increased temperature above 42 °C the ageing of a graft speeds up. Every vascular prosthesis is provided with production lot number and date of expiration.

DISPOSAL OF NON-USED MEDICAL DEVICE

Incineration together with municipal waste.

SIZE AND CATALOGUE CODE

Ra 1v K		
Length (mm)	Diameter (mm)	Catalogue code
300	6	022-06/300
	7	022-07/300
	8	022-08/300
	9	022-09/300
	10	022-10/300
	12	022-12/300
	14	022-14/300
	16	022-16/300
	18	022-18/300
	20	022-20/300
	22	022-22/300
	24	022-24/300
500	25	022-25/300
	26	022-26/300
	28	022-28/300
	30	022-30/300
	32	022-32/300
	6	022-0062
	7	022-0072
	8	022-0082

Ra 1v K		
Length (mm)	Diameter (mm)	Catalogue code
500	6	022-0062
	7	022-0072
	8	022-0082
	9	022-0092
	10	022-0102
	12	022-0122
	14	022-0142
	16	022-0162
	18	022-0182
	20	022-0202
	22	022-0222
	24	022-0242
	25	022-0252
	26	022-0262
	28	022-0282
	30	022-0302
	32	022-0322
	6	022-0062

Ra bv K		
Length (mm)	Diameter (mm)	Catalogue code
400	6/12	024-0612
	7/14	024-0714
	8/16	024-0816
	9/18	024-0918
	10/20	024-1020
	11/22	024-1122

SIZE MADE TO ORDER

Ra 1v K, Ra bv K	
Length 30–700 (mm)	Diameter 6–34 (mm)

EXPLANATORY NOTES ON THE SYMBOLS ON LABELS:		
Use by	Keep away from sunlight	
Catalogue code	Temperature limitation	
Sterilized by irradiation	Do not use if package is damaged	
Lot number (batch)	Waste bin	
Attention, see instructions for use!	Plastics	
Date of manufacture	Polyethyleneterephthalate / polyethylene	
Do not reuse	Conformity marking	
Producer	Do not re-sterilize	

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