Instructions for use

Pg	Language		Instructions for use
7	English	English	Instructions for use
9	French	Français	Mode d'emploi
11	German	Deutsch	Gebrauchsanweisung
13	Italian	Italiano	Istruzioni per l'uso
15	Spanish	Español	Instrucciones de uso
17	Dutch	Nederlands	Gebruiksaanwijzing
19	Portuguese	Português	Instruções de utilização
21	Greek	Ελληνικά	Οδηγίες χρήσης
24	Danish	Dansk	Brugsanvisning
26	Swedish	Svenska	Bruksanvisning
28	Finnish	Suomi	Käyttöohjeet
30	Norwegian	Norsk	Bruksanvisning
32	Polish	Polski	Instrukcja użytkowania
35	Hungarian	Magyar	Használati útmutató
37	Czech	Čeština	Návod k použití
39	Turkish	Türkçe	Kullanma talimatı
41	Russian	Русский	Инструкция по применению
44	Korean	한국어	사용법
46	Bulgarian	Български	Инструкции за употреба
48	Estonian	Eesti	Kasutusjuhend
50	Latvian	Latviski	Lietošanas pamācība
52	Lithuanian	Lietuviškai	Naudojimo instrukcija
54	Romanian	Română	Instrucțiuni de utilizare
56	Slovak	Slovenčina	Návod na použitie
58	Serbian	Srpski	Uputstvo za upotrebu

REF

Catalog Number | Numéro de catalogue | Bestellnummer | Numero di catalogo | Número de referencia | Catalogusnummer | Número de catálogo | Αριθμός καταλόγου | Katalognummer | Katalognummer | Luettelonumero | Katalognummer | Numer katalogowy | Katalógusszám | Katalogové číslo | Katalog Numarasi | Номер по каталогу | 카탈로그 번호 | Каталоген номер | Katalooginumber | Kataloga numurs | Katalogo numeris | Număr de catalog | Katalógové číslo | Kataloški broj

LOT

Batch Code | Code du lot | Chargenbezeichnung | Codice del lotto | Código de lote | Partijcode | Código de lote | Κωδικός παρτίδας | Batchkode | Batchkod | Eräkoodi | Partikode | Kod partii | Gyártási tétel kódja | Kód dávky | Parti Kodu | Код партии | 배치 코드 | Партиден код | Partii kood | Partijas kods | Partijos koda | Număr de lot | Kód šarže | Serijski broj

Balloon outer diameter | Diamètre externe du ballonnet | Ballonaufendurchmesser | Diametro esterno del palloncino | Diámetro externo del balón | Buitendiameter ballon | Diámetro externo do balão | Εξωτερική διάμετρος μπαλονιού | Udvendig ballondiameter | Ballongens ytterdiameter | Pallon ulkoläpimitta | Ytre ballongdiameter | Średnica zewnętrzna balonu | Külső ballonátmérő | Vnější průměr balónku | Balon dış çapı | Наружный диаметр баллона | 풍선 외경 | Външен диаметър на балона | Ballooni välisläbimõõt | Balona ārējais diametrs | Išorinis balionėlio skersmuo | Diametrul exterior al balonului | Vonkajši priemer balónika | Spoljni prečnik balona

Balloon Length | Longueur du ballonnet | Ballonlänge | Lunghezza del palloncino | Longitud del balón | Lengte ballon | Comprimento do balão | Μήκος μπαλονιού | Ballonlængde | Ballonglängd | Pallon pituus | Ballonglengde | Długość balonu | Ballonhossz | Délka balónku | Balon Uzunluğu | Длина баллона | 풍선 길이 | Дължина на балона | Ballooni pikkus | Balona garums | Balionėlio ilgis | Lungimea balonului | Dĺžka balónika | Dužina balona

NOM

Nominal Pressure | Pression nominale | Nenndruck | Pressione nominale | Presión nominal | Nominale druk | Pressão nominal | Ονομαστική πίεση | Nominelt tryk | Nominellt tryck | Nimellispaine | Nominelt trykk | Ciśnienie nominalne | Névleges nyomás | Jmenovitý tlak | Nominal Basınç | Номинальное давление | 표준 압력 | Номинално налягане | Nimirõhk | Nominālais spiediens | Vardinis slėgis | Presiune nominală | Menovitý tlak | Nominalni pritisak

RBP

Rated burst pressure | Pression de rupture | Nennberstdruck | Pressione nominale di scoppio | Presión nominal de rotura | Nominale barstdruk | Pressão estimada de rotura | Ονομαστική πίεση ρήξης | Nominelt bristetryk | Nominellt rupturtryck | Nimellinen puhkeamispaine | Nominelt sprengtrykk | Znamionowe ciśnienie rozerwania | Névleges repedési nyomás | Maximālnī poruchový tlak | Azami patlama basıncı (RBP) | Номинальное давление разрыва | 정격 파열 압력 | Номинално налягане на спукване | Nominaalne lõhkemisrõhk | Nominālais pārplīšanas spiediens | Vardinis plyšimo slėgis | Presiune de rupere nominală | Menovitý tlak prasknutia | Nominalni pritisak pucanja

Max guidewire compatibility | Compatibilité max. avec les guides métalliques | Max. Führungsdraht-kompatibilität | Compatibilità massima del filo guida | Compatibilidad máx. del alambre guía | Max. voerdraadcompatibiliteit | Compatibilidade com fio-guia | Συμβατότητα μέγ. οδηγού σύρματος | Maks. guidewirekompatibilitet | Maxstorlek på ledare för kompatibilitet | Ohjainlangan suurin yhteensopivuus | Maks. ledevaierkompatibilitet | Zgodność przewodnika (maks.) | Vezetődrót maximális kompatibilitása | Maximālnī kompatibilitā vodicīho drātu | Maks. klavuz tel uyumluluğu | Максимальный совместимый проводник | 최대 가이드와이어 호환성 | Макс. съвместимост на водач | Max juhtetraadiga ühilduvus | Maks. vadítástígas saderība | Maksimalus kreipiamosios vielos suderinamumas | Compatibilitatea maximă a firului de ghidaj | Maximālna kompatibilitā vodiaceho drātu | Maks. kompatibilnost žice vodiče

Min. guide catheter (MGC) | Cathéter-guide, min. | Min. Führungskatheter (MGC) | Diametro interno minimo del catetere guida (MGC) | Catéter guía mín. (MGC) | Min. geleidekatheter (MGC) | Cateter-guia mín. (MGC) | Ελάχιστη οδηγός καθετήρας (MGC) | Min. guidekatheter (MGC) | Minsta guidekatheter (MGC) | Pienin ohjainkatetri | Min. ledekatheter (MGC) | Min. wymiar cewnika prowadzącego | Vezetőkatéter legkisebb belső átmérője | Minimālnī vodicī katetr (MGC) | Min. klavuz kateter (MGC) | Минимальный проводниковый катетер (MGC) | 최소 가이드 카테터(MGC) | Мин. водещ катетър (MGC) | Min juhtekateeter (MGC) | Min. vadítőkathetrs (MGC) | Min. kreipiamasis kateteris (MGC) | Diametrul min. al cateterului de ghidaj (MGC) | Minimālny vodiaci katéter (MGC) | Min. veličina vodećeg katetera (MGC)

Use by date | Date de péremption | Verfallsdatum | Data di scadenza | Fecha de caducidad | Isterste gebruiksdatum | Prazo de validade | Ημερομηνία λήξης | Udløbsdato | Utgångsdatum | Käytettävä viimeistään | Brukes innen | Data przydatności do użycia | Szavatossági idő | Použit do data | Son kullanna tarihi | Годен до | 사용 기한 | Да се използва преди | Kasutamise lõpptähtpäev | Izlietot līdz | Sunaudoti iki nurodytos datos | Data expirării | Dátum expirácie | Rok upotrebe

Non-pyrogenic | Apyrogène | Nicht pyrogen |
 Non pirogeno | Apirógeno | Niet-pyrogeen |
 Não pirogênico | Μη πυρογόνος | Ikke-pyrogen |
 Ikke-pyrogen | Ei-pyrogeeninen | Ikke-pyrogen |
 Niepirogenny | Nem pirogén | Apyrogenní | Pirojenik
 değıldir | Апи́рогенно | 비발열성 | Непи́рогенно |
 Mittepirogeenne | Nepirogēns | Nepirogeninis |
 Nepirogen | Nepyrogénne | Nepirogeno

Upper limit of temperature | Limite supérieure
 de température | Obere Temperaturgrenze |
 Limite superiore di temperatura | Limite superior
 de temperatura | Bovenlimiet temperatuur |
 Limite superior da temperatura | Ανώτατο όριο
 θερμοκρασίας | Øvre temperaturgrænse | Øvre
 temperaturgräns | Lämpötilan yläraja | Øvre
 temperaturgrense | Górna granica temperatury |
 Felső hőmérsékleti határérték | Horní teplotní limit |
 Üst sıcaklık sınırı | Максимальная допустимая
 температура | 온도 상한치 | Горна граница за
 температура | Temperatuuri ülempiir | Temperatūras
 augšējā robeža | Viršutinė temperatūros riba | Limita
 superioară de temperatură | Horný teplotný limit |
 Gornja granica temperature

Importer | Importateur | Importeur | Importatore |
 Entidad importadora | Importeur | Importador |
 Εισαγωγέας | Importør | Importör | Maahantuojat |
 Importør | Importer | Importör | Dovoze | Íthalatc |
 Импортёр | 수입업자 | Вноситель | Importija |
 Importētājs | Importuotojas | Importator | Dovoza |
 Uvoznik

Packaging unit | Nombre d'unités par emballage |
 Verpackungseinheit | Pezzi per confezione |
 Unidades en envase | Verpakkingseenheid |
 Unidades por embalagem | Μονάδα συσκευασίας |
 Pakkeenhed | Förpackningsenhet | Pakkausyksikkö |
 Emballasjeenhet | Jednostka pakowania | Csomagolási
 egység | Jednotka balení | Paketleme ünitesi |
 Количество единиц в упаковке | 포장단위 | Единица
 упаковки | Pakendühik | Iesaiņojuma vienība |
 Vienetų pakuotėje | Unități per ambalaj | Jednotka
 balenia | Jedinica pakovanja

Keep dry | Tenir au sec | Vor Feuchtigkeit schützen |
 Conservare in luogo asciutto | Conservar en lugar
 seco | Droog bewaren | Manter seco | Να διατηρείται
 στεγνό | Opbevares tørt | Förvaras torrt | Säilytettävä
 kuivana | Må holdes tørt | Chronić przed wilgocią |
 Szárazon tartandó | Chránit před vlhkem | Kuru tutun |
 Хранить в сухом месте | 항상 건조한 상태로
 보관하십시오 | Пазете сухо | Hoida kuivana | Glabāt
 sausā vietā | Laikyti sausi | A se păstra la loc uscat |
 Uchovávať v suchu | Držati suvim

Pressure (atm/kPa) | Pression atmosphérique (atm/kPa) |
 Atmosphärendruck (atm/kPa) | Pressione atmosferica
 (atm/kPa) | Presión (atm/kPa) | Druk (atm/kPa) |
 Pressão (atm/kPa) | Πίεση (atm/kPa) | Tryk (atm/kPa) |
 Tryck atm (kPa) | Paine (atm/kPa) | Trykk (atm/kPa) |
 Ciśnienie (atm/kPa) | Nyomás (atm/kPa) |
 Tlak (atm/kPa) | Basinc (atm/kPa) | Давление (атм/кПа) |
 압력 (atm/kPa) | Hálaryenne (atm/kPa) | Rõhk (atm/kPa) |
 Spiediens (atm./kPa) | Ślęgis (atm/kPa) |
 Presiune (atm/kPa) | Tlak (atm/kPa) |
 Pritisak (u atmosferama/kPa)

Sterilized using ethylene oxide | Stérilisé à l'oxyde
 d'éthylène | Mit Ethylenoxid sterilisiert | Sterilizzato
 con ossido di etilene | Esterilizado con óxido de
 etileno | Gesteriliseerd met ethylenoxide | Esterilizado
 por óxido de etileno | Αποστειρωμένο με οξείδιο
 αιθυλενίου | Steriliseret med ethylenoxid | Steriliserad
 med etylenoxid | Steriloitu etyleenioksidilla |
 Sterilisert med etylenoksid | Wyjałowione tlenkiem
 etylenu | Etilén-oxidál sterilizálva | Sterilizováno
 ethylenoxidem | Etilen oksit kullanılarak sterilize
 edilmiştir | Стерилизовано оксидом этилена |
 산화 에틸렌 가스로 멸균 처리 | Стерилизовано
 с этиленов оксид | Steriliseeritud etüleenoksidiiga |
 Sterilizēts ar etilēnoksidu | Sterilizuota etileno
 oksidu | Sterilizat cu oxid de etilenă | Sterilizované
 etylénoksidom | Sterilizovano pomoću etilen oksida

Manufacturer | Fabricant | Hersteller | Produttore |
 Fabricante | Fabrikant | Fabricante | Κατασκευαστής |
 Producent | Tillverkare | Valmistaja | Produzent |
 Producent | Gyártó | Výrobce | Üretici |
 Производител | 제조사 | Производител | Tootja |
 Ražotājs | Gamintojas | Producător | Výrobca |
 Proizvođač

Shaft length | Longueur de la tige | Schaftlänge |
 Lunghezza dello stelo | Longitud del eje | Lengte
 schacht | Comprimento da haste | Μήκος άξονα |
 Skaftlængde | Skaftlängd | Varren pituus |
 Skaftlengde | Długość trzonu | Szár hossza | Délka
 tubusu | Şaft uzunluğu | Длина шфта | 축 길이 |
 Дължина на шфта | Varre pikkus | Kanála garums |
 Vamzdėlio ilgis | Lungimea tijei | Dĺžka drieku |
 Dužina osovine

Do not use if package is damaged | Ne pas utiliser si
 l'emballage est endommagé | Bei beschädigter Packung
 nicht verwenden | Non usare se la confezione è
 danneggiata | No utilizar el producto si el envase está
 dañado | Niet gebruiken als de verpakking beschadigd
 is | Não utilizar se a embalagem se encontrar
 danificada | Μη χρησιμοποιείτε το προϊόν εάν η
 συσκευασία έχει υποστεί ζημιά | Må ikke anvendes,
 hvis emballagen er beskadiget | Får ej användas om
 förpackningen är skadad | Ei saa käyttää, jos pakkaus on
 vaurioitunut | Må ikke brukes hvis pakken er skadet |
 Nie używać, jeśli opakowanie jest uszkodzone | Ne
 használja, ha a csomagolás sérült | Nepoužívat, jestliže
 je balení poškozeno | Ambalaj hasarlıysa kullanmayın |
 Не использовать, если упаковка повреждена |
 포장에 손상된 것은 사용하지 마십시오 | Да
 не се използва, ако опаковката е повредена |
 Mitte kasutada, kui pakend on kahjustunud |
 Nelietot, ja iesaiņojums ir bojāts | Nenaudoti, jei
 pakuotė pažeista | A nu se folosi în caz de deteriorare
 a ambalajului | Nepoužívajte, ak je obal poškodený |
 Nemojte koristiti ako je pakovanje oštećeno

Medical Device | Dispositif médical | Medizinprodukt |
 Dispositivo medico | Producto sanitario | Medisch
 hulpmiddel | Dispositivo médico | Ιατροτεχνολογικό
 προϊόν | Medicinsk udstyr | Medicinteknisk produkt |
 Lääkinnällinen laite | Medicinsk utstyr | Wyrób
 medyczny | Orvostechnikai eszköz | Zdravotnícký
 prostředek | Tibbi Cihaz | Устройство медицинского
 назначения | 의료 기기 | Медицинско изделие |
 Meditsiiniseade | Mediciniska ierice | Medicinos
 prietaisai | Dispozitiv medical | Zdravotnicka
 pomůcka | Medicinsko sredstvo

Single sterile barrier system | Système de barrière stérile unique | Einfaches Sterilbarrieresystem | Sistema a única barriera sterile | Sistema de barrera estéril única | Een enkel steriel barrièresysteem | Sistema de barreira estéril única | Σύστημα μονού στείρου φραγμού | Enkelt sterilt barrièresystem | Enkelt sterilbarriärsystem | Yksittäinen steriili estejärjestelmä | Enkelt sterilt barrieresystem | System pojedynczej bariery jałowej | Egyszeres sterilgátrendszer | Systém jednoduché sterilní bariéry | Tek steril bariyer sistemi | Однobarьерная система для защиты стерильности | 일회용 멸균장벽시스템 | Sistema с единична стерилна бариера | Ühekordne steriilne barjäärsüsteem | Viena sterilitātes aizsargslāņa sistēma | Vieno sterilaus barjero sistema | Sistem de barieră sterilă unică | Systém jednej sterilnej bariéry | Sistem pojedinačne sterilne barijere

Do not reuse | Ne pas réutiliser | Nicht wiederverwenden | Non riutilizzare | No reutilizar | Niet opnieuw gebruiken | Não reutilizar | Μη επαναχρησιμοποιείτε | Må ikke genbruges | Får ej återanvändas | Ei saa käyttää uudelleen | Må ikke brukes på nytt | Nie używać ponownie | Ne használja fel újra! | Nepoužívat opětvně | Tekrar kullanmayın | Не использовать, если упаковка повреждена | 재사용하지 마십시오 | Не употребляйте повторно | Mitte korduskasutada | Nelietot atkārtoti | Nenaudoti pakartotinai | A nu se reutiliza | Nepoužívajte opakovane | Ne koristiti ponovo

Manufactured for | Fabriqué pour | Hergestellt für | Prodotto per | Fabricado para | Vervaardigd voor | Fabricado para | Κατασκευάζεται για | Fremstillet for | Tillverkas för | Valmistuttaja | Produzent for | Wyprodukowano dla | Gyártás megrendelője | Vyrobeno pro | Şunun için üretilmiştir | Произведено для | 다음을 위해 제조됨 | Jedinečně vyrobeno za | Toodetud ettevõttele | Pasūtītājs | Pagaminta užsakius | Fabricat pentru | Vyrobené pre | Proizvedeno za

Date of manufacture / country of manufacture | Date de fabrication/Pays de fabrication | Herstellungsdatum/Herstellungsland | Data di fabbricazione / paese di fabbricazione | Fecha de fabricación/país de fabricación | Productiedatum/land van productie | Data de fabrico/país de fabrico | Ημερομηνία κατασκευής/Χώρα κατασκευής | Fremstillingsdato/fremstillingsland | Tillverkningsdatum/tillverkningsland | Valmistuspäivä/valmistus- maa | Produktionsdato/produktionsland | Data produkcji / kraj produkcji | Gyártás ideje /helye | Datum výroby / země výroby | Üretim tarihi / üretim ülkesi | Дата производства / страна производства | 제조일자 / 제조국 | Дата на производство / държава на производство | Tootmiskuupäev/tootmisriik | Izgatavošanas datums / izgatavošanas valsts ražotājs | Pagaminimo data / šalis | Data fabricației/Țara de fabricație | Datum výroby/krajina výroby | Datum proizvodnje / država proizvodnje

Unique device identifier | Identifiant unique du dispositif | Unique Device Identifier (eindeutige Produktkennung) | Identificativo unico del dispositivo | Identificador único del dispositivo | Uniek identificatiekenmerk van het hulpmiddel | Identificação única do dispositivo | Μοναδικό αναγνωριστικό τεχνολογικού προϊόντος | Unik enhedsidentifikator | Unik produktidentifisering (UDI) | Yksilöllinen laitetunniste | Unik enhedsidentifikator | Unikálny identyfikator wyrobu | Egyedi eszközazonosító | Jedinečný identifikátor prostředku | Benzersiz cihaz tanımlayıcı | Унікальний ідентифікатор устро́йства | 의료기기 고유식별코드 | Унікален ідентифікатор на изделието | Unikaalne seadme ID | Unikāls ierīces identifikators | Unikalus prietaiso identifikatorius | Identificator unic al dispozitivului | Jedinečné identifikačné číslo pomôcky | Jedinstvena oznaka sredstva

Do not resterilize | Ne pas restériliser | Nicht erneut sterilisieren | Non risterilizzare | No reesterilizar | Niet opnieuw steriliseren | Não reesterilizar | Μη επαναποστερίωνετε | Må ikke resteriliseres | Får ej resteriliseras | Ei saa steriloida uudelleen | Må ikke resteriliseres | Nie resterylizować | Ne sterilizálja újra! | Neprovádět opětovnou sterilizaci | Tekrar sterilize etmeyin | Не подвергать повторной стерилизации | 재살균하지 마십시오 | Да не се стерилизира повторно | Mitte resterilisierida | Nesterilizēt atkārtoti | Nesterilizuoti pakartotinai | A nu se resteriliza | Nesterilizujte opakovane | Ne sterilisati ponovo

Caution: do not exceed rated burst pressure | Attention : ne pas dépasser la pression de rupture | Achtung: Nennberstdruck nicht überschreiten | Attenzione! Non superare la pressione nominale di scoppio | Precución: No supere el valor de presión nominal de rotura | Let op: Nominale barstdruk niet overschrijden | Atenção: não exceder a pressão estimada de rotura | Προσοχή: Μην υπερβαίνετε την ονομαστική πίεση ρήξης | Forsigtigt: Det nominelle bristetryk må ikke overstiges | Försiktighet! Överskrid inte det nominella bristningstrycket | Huomio: nimellistä puhkeamispainetta ei saa ylittää | Forsiktig: Ikke overskrid nominelt sprengtrykk | Przestroga: nie przekraczać znamionowego ciśnienia rozerwania | Figyelem! Ne lépje túl a névleges repedési nyomást | Upozornění: Nepřekročte maximální poruchový tlak | Dikkat: Azami patlama basıncını aşmayın | Предупреждение: не превышайте номинальное давление разрыва | 주의: 정격 파열 압력을 초과하지 않을 것 | Внимание: не превышаяйте номинального налягание на спукване | Ettevaatust! Ärge ületage nominaalset lõhkemisrõhku | Uzmanību! Nepārsniedziet nominālo pārlīšanas spiedienu | Dmesio! Neviršykite vardinio plyšimo slėgio | Atentie: nu depășiți presiunea de rupere nominală | Upozornenie: Neprekračujte menovitý tlak prasknutia | Pažnja: nemojte primenjivati pritisak viši od nominalnog pritiska pucanja

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English- en RX PTCA scoring balloon catheter with hydrophilic coating

INSTRUCTIONS FOR USE

AngioSculpt Evo Percutaneous Transluminal Coronary Angioplasty (PTCA) Scoring Balloon Catheter Rapid Exchange (RX) Delivery System

NOTE: These instructions apply to all balloon diameters and lengths.

STERILE: Sterilized with ethylene oxide gas. **Do not use if the package is open or damaged.**

CONTENTS: One (1) AngioSculpt Evo RX PTCA scoring balloon catheter.

STORAGE: Store in a dry, dark, cool place.

DEVICE NAME

The device name is the AngioSculpt Evo scoring balloon catheter; the generic name is Coronary angioplasty balloon catheter, cutting/scoring.

DEVICE DESCRIPTION

The AngioSculpt Evo scoring balloon catheter is a standard PTCA catheter with a scoring balloon near the distal tip. The balloon is designed to expand to a specified diameter and length at a specified pressure. The distal end of the catheter has a conventional nylon-blend balloon and a nitinol scoring element with three spiral struts that wrap around the balloon. The struts create focal concentrations of dilating force, which minimizes balloon slippage and assists in the luminal expansion of stenotic coronary arteries. Conventional radiopaque markers aid in positioning the balloon in the stenosis.

Figure 1 shows the distal section of the catheter with the scoring balloon. The proximal end of the balloon is connected to a balloon inflation channel.


Figure 1: Distal section of AngioSculpt Evo scoring balloon catheter

The catheter is compatible with standard 0.014" coronary guide wires and 5F guiding catheters. The catheter length is approximately 139 cm and is offered in diameters of 2.0–3.5 mm in 0.5 mm increments. The working length of the scoring element ranges from 6 to 20 mm. The catheter is supplied sterile and is intended for a single use. Figure 2 depicts the AngioSculpt catheter with the hydrophilic coating segment shown.

Figure 2: AngioSculpt Evo scoring balloon catheter with hydrophilic coating

INTENDED PURPOSE

The device is intended for balloon dilation of lesions in the coronary vasculature that are not easily dilatable with a standard angioplasty balloon.

INTENDED USERS

Intended users of the AngioSculpt Evo scoring balloon catheter are physicians trained to perform percutaneous transluminal coronary angioplasty and associated interventional procedures.

INTENDED PATIENT POPULATION

The intended patient population are those presenting with coronary artery disease.

INDICATIONS FOR USE

The AngioSculpt Evo scoring balloon catheter is indicated for the treatment of a hemodynamically significant coronary artery stenosis, including in-stent restenosis.

CLINICAL BENEFITS

The intended clinical benefit of the AngioSculpt Evo scoring balloon catheter is facilitation of complete stent expansion in lesions not amenable to balloon dilation prior to stent placement or to increase the residual lumen in ISR lesions in the coronary vasculature.

CONTRAINDICATIONS

The AngioSculpt Evo scoring balloon catheter should not be used for the following:

- Coronary artery lesions unsuitable for treatment by percutaneous revascularization.
- Coronary artery spasm in the absence of a significant stenosis.

WARNINGS

- The AngioSculpt Evo scoring balloon catheter is indicated for use in coronary arteries. The safety and effectiveness of the device has not been established, or is unknown, in vascular regions other than those specifically indicated.
- The AngioSculpt Evo scoring balloon catheter is coated with a hydrophilic coating on the transition tube at the distal end of the catheter for a length of 53 mm; see Figure 2. Please refer to Directions for Use, Step 16 for further information on how to prepare and use this catheter to ensure it performs as intended. Failure to abide by the warnings in this labeling might result in damage to the device coating, which may necessitate intervention.
- Administer appropriate antiplatelet, anticoagulant and coronary vasodilator therapy, consistent with institutional practice for coronary stent procedures, during and after the procedure.
- This device is intended for single (one) use only. Do not resterilize and/or reuse, as this can potentially result in compromised device performance and increased risk of inappropriate resterilization and cross contamination.
- For use in de novo or in-stent restenosis (ISR) lesions, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis. When used to pre-dilate the lesion prior to pre-planned stenting, the catheter should be one size smaller than the estimated vessel diameter (e.g., a 2.5 mm diameter device should be used in a vessel estimated to have a 3.0 mm diameter).
- PTCA in patients who are not acceptable candidates for coronary artery bypass graft surgery requires careful consideration, including possible hemodynamic support during PTCA, as treatment of this patient population carries special risk.
- When the catheter is exposed to the vascular system, it should be manipulated while under high quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.
- Balloon pressure should not exceed the rated burst pressure (RBP). [The RBP is based on results of in-vitro testing. At least 99.9% of the balloons (with 95% confidence) will not burst at or below their RBP. Use of a pressure monitoring device is recommended to prevent over-pressurization.]
- PTCA with the AngioSculpt Evo scoring balloon catheter should only be performed at hospitals where emergency coronary artery bypass graft surgery can be quickly performed on site (or at a nearby facility) in the event of a potentially injurious or life-threatening complication.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Proceed cautiously when using the AngioSculpt Evo scoring balloon catheter in a freshly deployed bare metal or drug eluting stent. The AngioSculpt catheter has not been tested for post-dilation of stents or lesions distal to freshly deployed stents in clinical studies. Bench testing has shown no additional risk when inserting or withdrawing the AngioSculpt catheter through stents (no interference with stent struts, no retention of or damage

to AngioSculpt catheter).

- Use the catheter prior to the "Use by date" specified on the package.

PRECAUTIONS

- Prior to angioplasty, the catheter should be examined to verify functionality, device integrity and to ensure that its size and length are suitable for the specific procedure for which it is to be used.
- Only physicians trained in the performance of percutaneous transluminal coronary angioplasty should use the catheter system.
- Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained.
- Do not rotate the catheter luer hub in excess of five (5) turns during use.
- Do not advance or retract the AngioSculpt Evo scoring balloon catheter over the floppy portion of the guide wire.
- Catheter manipulation, including advancement and retraction, should be performed by grasping the hypotube shaft.
- If unusual resistance is felt when the catheter is being manipulated or if it is suspected that the guide wire has become kinked, carefully remove the entire catheter system (AngioSculpt catheter and steerable guide wire) as a unit.
- If fluoroscopic guidance indicates that the AngioSculpt Evo scoring balloon catheter has advanced beyond the end of the guide wire, withdraw the catheter and reload the wire before advancing again.

ADVERSE EFFECTS

Possible complications include, but are not limited to, those listed in Table 1 below:

Table 1: Possible Complications

- Death
- Heart attack (acute myocardial infarction)
- Embolism
- Total occlusion of the treated coronary artery
- Coronary artery dissection, perforation, rupture, or injury
- Pericardial tamponade
- No/slow reflow of treated vessel
- Emergency coronary artery bypass (CABG)
- Emergency percutaneous coronary intervention
- CVA/stroke/embolic stroke
- Pseudoaneurysm
- Restenosis of the dilated vessel
- Unstable chest pain (angina)
- Thromboembolism or retained device components
- Irregular heart rhythm (arrhythmias, including life-threatening ventricular fibrillation)
- Severe low (hypotension)/high (hypertension) blood pressure
- Coronary artery spasm
- Hemorrhage or hematoma
- Need for blood transfusion
- Surgical repair of vascular access site
- Creation of a pathway for blood flow between the artery and the vein in the groin (arteriovenous fistula)
- Drug reactions, allergic reactions to x-ray dye (contrast medium)
- Infection

SPECIFICATION OF CLINICAL BENEFITS

The AngioSculpt Evo scoring balloon catheter summary of clinical studies/specification of clinical benefits is provided in the summary of safety and clinical performance (SSCP) which will be available in the European database on medical devices (EUDAMED) at the following URL: <https://ec.europa.eu/tools/eudamed>. The SSCP is linked to the GTIN, the 14-digit device identifier (UDI-DI) following the (01) within the UDI on the device label.

MATERIALS REQUIRED FOR USE WITH THE ANGIOSCULPT EVO CATHETER

WARNING - Use single use items only. Do not resterilize or reuse.

- Femoral, brachial or radial guiding catheter (≥ 5F)
- Hemostatic valve
- Contrast medium diluted 1:1 with normal saline
- Sterile heparinized normal saline
- 10 cc and 20 cc syringes for flushing and balloon prep
- Inflation device (indeflator)
- 0.014" coronary guide wire
- Guide wire introducer
- Guide wire torque device
- Radiographic contrast
- Manifold (for pressure monitoring and contrast injection), extension pressure tubing
- Sterile gauze

INSTRUCTIONS FOR USE

Prior to use of the AngioSculpt Evo Scoring balloon catheter, examine carefully for damage and device integrity. Do not use if the catheter has bends, kinks, missing components or other damage.

1. Pre-medicate patients with ASA, Clopidogrel/Ticlopidine, intravenous anti-coagulants, coronary vasodilators and GP2b/3a blockers according to institutional protocol for percutaneous coronary interventions involving stents.
2. Perform coronary angiogram in at least 2 orthogonal views of the target lesion prior to device deployment with online QCA.
3. Position 0.014" coronary guide wire of choice beyond the target lesion.
4. In the case of pre-planned stenting, select a catheter that is one size smaller relative to the estimated reference vessel diameter (RVD). For example, if the RVD by visual "estimation" during the procedure is 3.0 mm, use a 2.5 mm diameter AngioSculpt Evo scoring balloon catheter. For in-stent restenosis, select a catheter that closely approximates the RVD ($\leq 1.0 \times \text{RVD}$).
5. Using sterile technique, remove an appropriately sized ($\leq 1.0 \times$ reference vessel diameter (RVD)) AngioSculpt device from the sterile package and place on the sterile field.
6. Remove the stylet from the catheter.
7. Remove the protective tubing from the balloon.
8. Inspect the device to ensure that all components are intact.
9. Flush the guide wire lumen with saline by carefully inserting the distal catheter tip into the distal end of a 10 cc syringe and injecting saline until droplets emerge from the proximal guide wire lumen.
10. Attach a stopcock to the catheter's balloon inflation port.
11. Attach 20 cc syringe filled with 2-3 cc of radiographic contrast to the catheter balloon inflation port.
12. Open the stopcock to the syringe, aspirate/remove air from the catheter balloon lumen using the 20 cc syringe filled with 2-3 cc of radiographic contrast and leave on vacuum for 30 seconds.
13. Close the stopcock to the catheter balloon inflation port and remove the syringe.
14. Attach inflation device (indeflator), filled with 50:50 mixture of radiographic contrast and normal saline, to the balloon inflation port by creating a meniscus. Avoid introducing air bubbles into the catheter balloon lumen.
15. Open the stopcock to the inflation device and aspirate using the inflation device, locking in vacuum.
NOTE: All air must be removed from the balloon and displaced with contrast medium prior to inserting into the body (repeat steps 11-14, if necessary).
16. To activate the hydrophilic coating, wet the catheter's coated area with saline or wipe a wet sterile gauze over the coated area immediately prior to insertion into the guide catheter (see Figure 2).

NOTE: Avoid using alcohol, antiseptic solutions, or other solvents to wet the hydrophilic coating, because this may cause unpredictable changes in the coating which could affect the device safety and performance.

17. Advance the AngioSculpt Evo scoring balloon catheter over the coronary guide wire to the target lesion (may predilate with a 2.0 – 2.5 mm diameter balloon catheter if needed to cross the lesion). **NOTE: When backloading the catheter onto the guide wire, the catheter should be supported, ensuring that the guide wire does not come in contact with the balloon.**
18. Inflate the AngioSculpt Evo balloon per the following recommended protocol:
 - Increase the inflation pressure by 2 atmospheres every 10-15 seconds until full device inflation is achieved
 - May inflate to a maximum pressure that is $\leq$ RBP at the discretion of the physician (bearing in mind the estimated inflated diameter of the device at a given pressure)
19. Perform coronary angiogram (in the same view(s) as step 2) of the target lesion following completion of device treatment (and prior to adjunctive stent placement).
20. To remove the AngioSculpt Evo scoring balloon catheter, apply negative pressure to the inflation device and confirm that the balloon is fully deflated. The catheter should be retracted only by grasping the hypotube shaft. If there is difficulty removing the catheter through the hemostatic valve of the guide catheter due to lubricity, gently grasp the catheter with a wet sterile gauze. **NOTE: Do not rotate the catheter shaft in excess of 180 degrees when the tip is constrained. Do not rotate the catheter luer hub in excess of five (5) turns during use. Catheter manipulation, including advancement and retraction, should be performed by grasping the catheter shaft.**
21. Inspect all components to ensure that the device is intact. Follow institutional procedures for disposal of biohazards. If device malfunction occurs or any defects are noted on inspection, flush the guide wire lumen and clean the outer surface of the device with saline, store the device in a sealed plastic bag and contact Philips or the EC Authorized Representative for further instructions.
22. Complete any additional interventions as clinically indicated (e.g. stent placement).
23. Remove the coronary guide wire and perform coronary angiography (in the same view(s) as step 2) of the target lesion following completion of all interventions.
24. Remove all catheters and manage the arterial access site according to institutional protocol.
25. Continue treatment with ASA, Clopidogrel/Ticlopidine, and GP2b/3a blockers according to institutional protocol for percutaneous coronary interventions involving stents

REPORTING OF A SERIOUS INCIDENT, CUSTOMER COMPLAINTS AND FEEDBACK

If a serious incident (death or any intervention) has occurred and the AngioSculpt Evo scoring balloon catheter was in use at any time, it should be reported to the manufacturer immediately and/or the competent authority of the Member state in which the user and/or patient is established. A serious incident means any incident that directly or indirectly led, might have led or, in case of recurrence, could lead to any of the following:

- Intervention to prevent a serious injury
- The death of a patient, user or other person
- The temporary or permanent serious deterioration of a patient's, user's, fetus or other person's state of health
- A serious public health threat

Any complaints against the devices that alleges deficiencies

related to the identity, quality, durability, reliability, safety, effectiveness, performance, or conformance to a specification of a Philips medical device, should be reported to the post market department at SPNC-Complaints@philips.com.

Any other device feedback can be reported to the post market department at SPNC-Complaints@philips.com.

REFERENCES

The physician should consult recent literature on current medical practice regarding balloon dilatation and PTCA procedures.

DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

There is no express or implied warranty, including any implied warranty of merchantability or fitness for a particular purpose, on the manufacturer product(s) described in this publication. Under no circumstances shall the manufacturer be liable for any direct, incidental or consequential damages other than as expressly provided by specific law. No person has the authority to bind the manufacturer to any representation or warranty except as specifically set forth herein.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

The AngioSculpt Evo scoring balloon catheter summary of safety and clinical performance (SSCP) will be available in the European database on medical devices (EUDAMED) at the following URL: <https://ec.europa.eu/tools/eudamed>. The SSCP is linked to the GTIN, the 14-digit device identifier (UDI-DI) following the (01) within the UDI on the device label.

DISPOSAL

After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

REACH DECLARATION

REACH requires Philips Healthcare (PH) to provide chemical content information for "substances of Very High Concern (SVHC)" if they are present above 0.1% of the product weight." The SVHC is updated on a regular basis. Therefore, refer to the following Philips REACH website for the most up-to-date list of products containing SVHC above the threshold: <http://www.philips.com/about/sustainability/reach.page>

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Philips Medical Systems Nederland B.V.
Philips Electronics Australia Ltd
Philips Image Guided Therapy Corporation